

Product datasheet for **HP208411**

ETS1 Human qPCR Primer Pair (NM_005238)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2113
Forward Sequence:	GAGTCAACCCAGCCTATCCAGA
Reverse Sequence:	GAGCGTCTGATAGGACTCTGTG
ACCN:	NM_005238 , NM_005238.1 , NM_005238.2 , NM_005238.3 , BC017314 , BG035366 , BM468663 , BP365215 , BT019452 , BT019453 , BX640634 , NM_005238.4
UniProt ID:	P14921
Synonyms:	c-ets-1; ETS-1; EWSR2; p54
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)