

Product datasheet for **HP208424**

FRA2 (FOSL2) Human qPCR Primer Pair (NM_005253)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2355
Forward Sequence:	AAGAGGAGGAGAAGCGTCGCAT
Reverse Sequence:	GCTCAGCAATCTCCTTCTGCAG
ACCN:	NM_005253 , NM_005253.1 , NM_005253.2 , NM_005253.3 , BC022791 , BC022791.1 , BC008899 , BC085197 , BM023443 , BQ574214 , BX647822 , NM_005253.4
UniProt ID:	P15408
Synonyms:	FRA2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)