

## Product datasheet for **HP208430**

### GDF 9 (GDF9) Human qPCR Primer Pair (NM\_005260)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2661
Forward Sequence:	GTCACCTCTACAACACTGTTCGG
Reverse Sequence:	CAACGGTAGTAATGCGATCCAGG
ACCN:	<a href="#">NM_005260</a> , <a href="#">NM_005260.1</a> , <a href="#">NM_005260.2</a> , <a href="#">NM_005260.3</a> , <a href="#">NM_005260.4</a> , <a href="#">NM_005260.5</a> , <a href="#">BC096228</a> , <a href="#">BC096229</a> , <a href="#">BC096230</a> , <a href="#">BC096231</a> , <a href="#">BX117145</a>
UniProt ID:	<a href="#">O60383</a>
Synonyms:	POF14
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)