

Product datasheet for **HP208444**

GNL1 Human qPCR Primer Pair (NM_005275)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2794
Forward Sequence:	GCAGACTGATTCAGCAATGGAGC
Reverse Sequence:	GGTTCTGGAGACACTCAGACT
ACCN:	NM_005275 , NM_005275.1 , NM_005275.2 , NM_005275.3 , BC018366 , BC018366.1 , BC006343 , BC013959 , BC063657 , BM833876 , BQ188153 , BT006648 , BU687484
UniProt ID:	P36915
Synonyms:	HSR1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)