

## Product datasheet for **HP208452**

### GPR12 Human qPCR Primer Pair (NM\_005288)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2835
Forward Sequence:	GCTTCAGTCAGAAGCCACCAAG
Reverse Sequence:	TACAGTGAGAGGTAGCGGTCAAC
ACCN:	<a href="#">NM_005288</a> , <a href="#">NM_005288.1</a> , <a href="#">NM_005288.2</a> , <a href="#">NM_005288.3</a> , <a href="#">BC112144</a> , <a href="#">BC112144.1</a> , <a href="#">BC067448</a> , <a href="#">BC067449</a> , <a href="#">BC067450</a> , <a href="#">BC067451</a> , <a href="#">BC067452</a> , <a href="#">BC112146</a> , <a href="#">BM727022</a> , <a href="#">NM_005288.4</a>
UniProt ID:	<a href="#">P47775</a>
Synonyms:	GPCR12; GPCR21; PPP1R84
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online >](#)