

Product datasheet for **HP208467**

GPR43 (FFAR2) Human qPCR Primer Pair (NM_005306)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2867
Forward Sequence:	GCCTGGGTTATGTCCTTTGGTC
Reverse Sequence:	CGGTGAAGTTCTCGTAGCAGGT
ACCN:	<u>NM_005306</u> , <u>NM_005306.1</u> , <u>NM_005306.2</u> , <u>BC095535</u> , <u>BC096198</u> , <u>BC096199</u> , <u>BC096200</u> , <u>BC096201</u> , <u>NM_005306.3</u>
UniProt ID:	<u>O15552</u>
Synonyms:	FFA2R; GPR43
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)