

Product datasheet for **HP208483**

Histone H1.1 (HIST1H1A) Human qPCR Primer Pair (NM_005325)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3024
Forward Sequence:	GGCAAAGAAACCTGCTAAGGCTG
Reverse Sequence:	TAAGAGCTGCCAACGACACACC
ACCN:	<u>NM_005325</u> , <u>NM_005325.1</u> , <u>NM_005325.2</u> , <u>NM_005325.3</u> , <u>BC112140</u> , <u>BC112140.1</u> , <u>BC069492</u> , <u>BC101593</u> , <u>NM_005325.4</u>
UniProt ID:	<u>Q02539</u>
Synonyms:	H1.1; H1A; H1F1; HIST1; HIST1H1A
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)