

## Product datasheet for **HP208491**

### HCLS1 Human qPCR Primer Pair (NM\_005335)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3059
Forward Sequence:	GCTATGGAAGAGCCAGCAGTAC
Reverse Sequence:	TGCTGGTTCTCACAGGCTCAGA
ACCN:	<a href="#">NM_005335</a> , <a href="#">NM_005335.1</a> , <a href="#">NM_005335.2</a> , <a href="#">NM_005335.3</a> , <a href="#">NM_005335.4</a> , <a href="#">NM_005335.5</a> , <a href="#">BC016758</a> , <a href="#">BC016758.2</a> , <a href="#">BP338037</a> , <a href="#">BQ002360</a> , <a href="#">BT006824</a> , <a href="#">BY994733</a>
UniProt ID:	<a href="#">P14317</a>
Synonyms:	CTTNL; HS1; lckBP1; p75
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)