

## Product datasheet for **HP208502**

### **RBPJK (RBPJ) Human qPCR Primer Pair (NM\_005349)**

#### **Product data:**

<b>Product Type:</b>	qPCR Primer Pairs
<b>Gene ID:</b>	3516
<b>Forward Sequence:</b>	TCATGCCAGTTCACAGCAGTGG
<b>Reverse Sequence:</b>	TGGATGTAGCCATCTCGGACTG
<b>ACCN:</b>	<a href="#">NM_005349</a> , <a href="#">NM_005349.1</a> , <a href="#">NM_005349.2</a> , <a href="#">NM_005349.3</a> , <a href="#">BC020780</a> , <a href="#">BC064976</a> , <a href="#">BC064976.1</a> , <a href="#">BC053531</a> , <a href="#">BC053568</a> , <a href="#">BE269143</a> , <a href="#">BM472788</a> , <a href="#">BM555733</a> , <a href="#">BM560796</a> , <a href="#">BM680668</a> , <a href="#">BQ216033</a> , <a href="#">NM_005349.4</a>
<b>UniProt ID:</b>	<a href="#">Q06330</a>
<b>Synonyms:</b>	AOS3; CBF-1; CBF1; csl; IGKJRB; IGKJRB1; KBF2; RBP-J; RBP-JK; RBP-J kappa; RBPJK; RBPSUH; SUH
<b>Components:</b>	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
<b>Quality Control:</b>	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
<b>Storage:</b>	Store at -20°C.
<b>Stability:</b>	The primer mix is stable for one year from date of shipping.



[View online »](#)