

Product datasheet for **HP208515**

Myoglobin (MB) Human qPCR Primer Pair (NM_005368)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4151
Forward Sequence:	AGTCGCATGCCACCAAGCACAA
Reverse Sequence:	CCTGGGCATCAGCACCAAAGTC
ACCN:	NM_005368 , NM_005368.1 , NM_005368.2 , BC014547 , BC014547.1 , BC018001 , BQ925917 , BQ956082 , BU585249 , NM_005368.3
UniProt ID:	P02144
Synonyms:	PVALB
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)