

Product datasheet for **HP208517**

METTL1 Human qPCR Primer Pair (NM_005371)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4234
Forward Sequence:	TGGCTTCCAGAACATCGCCTGT
Reverse Sequence:	TGTCCGCTTGAAATGTGGGTCG
ACCN:	<u>NM_005371</u> , <u>NM_005371.1</u> , <u>NM_005371.2</u> , <u>NM_005371.4</u> , <u>NM_005371.5</u> , <u>BC000550</u> , <u>BG420288</u> , <u>BQ651940</u> , <u>NM_005371.6</u>
UniProt ID:	<u>Q9UBP6</u>
Synonyms:	C12orf1; TRM8; TRMT8; YDL201w
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)