

Product datasheet for **HP208529**

Neuronatin (NNAT) Human qPCR Primer Pair (NM_005386)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4826
Forward Sequence:	GTGTTCCCTGGAATGCTGCATTTAC
Reverse Sequence:	GACACCGTGTATGCCAGCTTCT
ACCN:	NM_005386 , NM_005386.1 , NM_005386.2 , NM_005386.3 , BC001768 , BC001768.1 , BE255132 , BG472075 , BT019426 , NM_005386.4
UniProt ID:	Q16517
Synonyms:	Peg5
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)