

Product datasheet for **HP208539**

AMPK beta 2 (PRKAB2) Human qPCR Primer Pair (NM_005399)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5565
Forward Sequence:	GAAGGAGCACAAGATCATGGTGG
Reverse Sequence:	TACGGAGTCCTCCAAATCCTGC
ACCN:	<u>NM_005399</u> , <u>NM_005399.1</u> , <u>NM_005399.2</u> , <u>NM_005399.3</u> , <u>BC053610</u> , <u>BC053610.1</u> , <u>BM563523</u> , <u>BM724633</u> , <u>BM769516</u> , <u>BU166852</u> , <u>BU619639</u> , <u>BX478507</u> , <u>NM_005399.5</u>
UniProt ID:	<u>O43741</u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)