

Product datasheet for **HP208550**

SnoN (SKIL) Human qPCR Primer Pair (NM_005414)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6498
Forward Sequence:	GCAGGAAGGTGACCATGTTTCTC
Reverse Sequence:	TACTTGCCTCTGTCTTTGTGAGC
ACCN:	NM_005414 , NM_005414.1 , NM_005414.2 , NM_005414.3 , NM_005414.4 , BC059386 , BC059386.1 , BC041913
UniProt ID:	P12757
Synonyms:	SNO; SnoA; SnoI; SnoN
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)