

Product datasheet for **HP208554**

STAT2 Human qPCR Primer Pair (NM_005419)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6773
Forward Sequence:	CAGGTCACAGAGTTGCTACAGC
Reverse Sequence:	CGGTGAACTTGCTGCCAGTCTT
ACCN:	<u>NM_005419</u> , <u>NM_005419.1</u> , <u>NM_005419.2</u> , <u>NM_005419.3</u> , <u>BC051284</u> , <u>BC051284.1</u> , <u>BC006092</u> , <u>BX640607</u> , <u>NM_005419.4</u>
UniProt ID:	<u>P52630</u>
Synonyms:	IMD44; ISGF-3; P113; PTORCH3; STAT113
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)