

Product datasheet for **HP208568**

CCDC6 Human qPCR Primer Pair (NM_005436)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8030
Forward Sequence:	GAACTTGAGGCTCCAGAGGAAG
Reverse Sequence:	CTGGACACAGTGCGAGGTCTTA
ACCN:	NM_005436 , NM_005436.1 , NM_005436.2 , NM_005436.3 , NM_005436.4 , BC036757 , BC036757.2 , BC022284 , BC064391 , BM672399 , NM_005436.5
UniProt ID:	Q16204
Synonyms:	D10S170; H4; PTC; TPC; TST1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)