

Product datasheet for **HP208597**

NAALAD2 Human qPCR Primer Pair (NM_005467)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10003
Forward Sequence:	GGGCATCAACTGTACTGGGAAG
Reverse Sequence:	GGATATGGCTGTACCTCAGGAG
ACCN:	<u>NM_005467</u> , <u>NM_005467.2</u> , <u>NM_005467.3</u> , <u>BC096316</u> , <u>BC027594</u> , <u>BC038840</u> , <u>BC067793</u> , <u>BC096317</u> , <u>BC096318</u> , <u>BC099646</u> , <u>BU618706</u> , <u>NM_005467.4</u>
UniProt ID:	<u>Q9Y3Q0</u>
Synonyms:	GPCIII; GPCIII; NAADALASE2; NAALADASE2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)