

Product datasheet for **HP208603**

GLCNE (GNE) Human qPCR Primer Pair (NM_005476)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10020
Forward Sequence:	GAAGCATACGCCTCTGGAATGG
Reverse Sequence:	CAGCCTCATCTTTGGCACTGAC
ACCN:	NM_005476 , NM_005476.1 , NM_005476.2 , NM_005476.3 , NM_005476.4 , NM_005476.5 , BC121179 , BF197948 , BM996750 , BM999446 , BP201585 , BQ217959 , NM_005476.7
UniProt ID:	Q9Y223
Synonyms:	DMRV; GLCNE; IBM2; NM; Uae1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)