

Product datasheet for **HP208622**

ABCA1 Human qPCR Primer Pair (NM_005502)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	19
Forward Sequence:	CAGGCTACTACCTGACCTTGGT
Reverse Sequence:	CTGCTCTGAGAAACACTGTCCTC
ACCN:	BC034824 , NM_005502 , NM_005502.1 , NM_005502.2 , NM_005502.3 , BC146856 , BU753834 , BX955241 , NM_005502.4
UniProt ID:	O95477
Synonyms:	ABC-1; ABC1; CERP; HDLCQTL13; HDLDT1; HPALP1; TGD
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)