

## Product datasheet for **HP208625**

### LIMPII (SCARB2) Human qPCR Primer Pair (NM\_005506)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	950
Forward Sequence:	GCCAATACGTCAGACAATGCCG
Reverse Sequence:	CTCATCTGCTTGGTAAAAGTGTGG
ACCN:	<u><a href="#">NM_005506</a></u> , <u><a href="#">NM_005506.1</a></u> , <u><a href="#">NM_005506.2</a></u> , <u><a href="#">NM_005506.3</a></u> , <u><a href="#">BC021892</a></u> , <u><a href="#">BC021892.1</a></u> , <u><a href="#">BC035111</a></u> , <u><a href="#">BQ924546</a></u> , <u><a href="#">BT006939</a></u> , <u><a href="#">BT007661</a></u> , <u><a href="#">NM_005506.4</a></u>
UniProt ID:	<u><a href="#">Q14108</a></u>
Synonyms:	AMRF; CD36L2; EPM4; HLGP85; LGP85; LIMP-2; LIMPII; SR-BII
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)