

Product datasheet for **HP208626**

Cofilin 1 (CFL1) Human qPCR Primer Pair (NM_005507)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1072
Forward Sequence:	GCAACCTATGAGACCAAGGAGAG
Reverse Sequence:	TCTTGATGGCGTCCTTGGAGCT
ACCN:	<u>NM_005507</u> , <u>NM_005507.1</u> , <u>NM_005507.2</u> , <u>BC012318</u> , <u>BC012318.1</u> , <u>BC011005</u> , <u>BC012265</u> , <u>BC018256</u> , <u>BC024975</u> , <u>BT006846</u> , <u>NM_005507.3</u>
UniProt ID:	<u>P23528</u>
Synonyms:	CFL; cofilin; HEL-S-15
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)