

Product datasheet for **HP208633**

HLAB (HLA-B) Human qPCR Primer Pair (NM_005514)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3106
Forward Sequence:	CTGCTGTGATGTGTAGGAGGAAG
Reverse Sequence:	GCTGTGAGAGACATCAGAGC
ACCN:	<u>BC013187</u> , <u>NM_005514</u> , <u>NM_005514.1</u> , <u>NM_005514.2</u> , <u>NM_005514.3</u> , <u>NM_005514.4</u> , <u>NM_005514.5</u> , <u>NM_005514.6</u> , <u>NM_005514.7</u> , <u>BC013187.1</u> , <u>BC091497</u> , <u>BC007243</u> , <u>BC013828</u> , <u>BX647315</u>
UniProt ID:	<u>P01889</u>
Synonyms:	AS; B-4901; HLAB
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)