

Product datasheet for **HP208648**

Heparan Sulfate Proteoglycan 2 (HSPG2) Human qPCR Primer Pair (NM_005529)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3339
Forward Sequence:	TCAGGCGAGTATGTGTGCCATG
Reverse Sequence:	GATGAAGACTCGATCCTGACAGG
ACCN:	NM_005529 , NM_005529.1 , NM_005529.2 , NM_005529.3 , NM_005529.4 , NM_005529.5 , NM_005529.6 , BC033152 , BC109204 , BE273742 , BG056875
UniProt ID:	P98160
Synonyms:	HSPG; PLC; PRCAN; SJA; SJS; SJS1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)