

Product datasheet for **HP208653**

IFNGR2 Human qPCR Primer Pair (NM_005534)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3460
Forward Sequence:	CTCCATTCTGCCTGGGTGACAA
Reverse Sequence:	CGTGGAGGTATCAGCGATGTCA
ACCN:	<u>NM_005534</u> , <u>NM_005534.1</u> , <u>NM_005534.2</u> , <u>NM_005534.3</u> , <u>BC003624</u> , <u>BG391103</u> , <u>BQ008968</u> , <u>NM_005534.4</u>
UniProt ID:	<u>P38484</u>
Synonyms:	AF-1; IFGR2; IFNGT1; IMD28
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)