

## Product datasheet for **HP208654**

### IL12RB1 Human qPCR Primer Pair (NM\_005535)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3594
Forward Sequence:	TGAGATTCTCGGTGGAGCAGCT
Reverse Sequence:	CTGTAGTCGGTAAGTGACCTCC
ACCN:	<u><a href="#">NM_005535</a></u> , <u><a href="#">NM_005535.1</a></u> , <u><a href="#">NM_005535.2</a></u> , <u><a href="#">BC029121</a></u> , <u><a href="#">BC137404</a></u> , <u><a href="#">BC137406</a></u> , <u><a href="#">BX647221</a></u> , <u><a href="#">NM_005535.3</a></u>
UniProt ID:	<u><a href="#">P42701</a></u>
Synonyms:	CD212; IL-12R-BETA1; IL12RB; IMD30
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)