

## Product datasheet for **HP208657**

### Inhibin beta C chain (INHBC) Human qPCR Primer Pair (NM\_005538)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3626
Forward Sequence:	CTGTGTCCAGAGCTGCTTTGAG
Reverse Sequence:	AGACGAGTCTGGTTGATGGTGG
ACCN:	<u><a href="#">NM_005538</a></u> , <u><a href="#">NM_005538.1</a></u> , <u><a href="#">NM_005538.2</a></u> , <u><a href="#">NM_005538.3</a></u> , <u><a href="#">BC020693</a></u> , <u><a href="#">BC130324</a></u> , <u><a href="#">BC130326</a></u> , <u><a href="#">BP237373</a></u> , <u><a href="#">BX105206</a></u> , <u><a href="#">NM_005538.4</a></u>
UniProt ID:	<u><a href="#">P55103</a></u>
Synonyms:	IHBC
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)