

Product datasheet for **HP208658**

INPP5A Human qPCR Primer Pair (NM_005539)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3632
Forward Sequence:	AGCGAGGAGAAGGTTGTACCT
Reverse Sequence:	CATGGGCATGAGGTTTACCTGC
ACCN:	NM_005539 , NM_005539.1 , NM_005539.2 , NM_005539.3 , NM_005539.4 , BC096280 , BC052265 , BC062300 , BC096281 , BC117515 , BI964266 , BM467443 , BM666428 , NM_005539.5
UniProt ID:	Q14642
Synonyms:	5PTASE
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)