

Product datasheet for **HP208670**

KRTAP5-9 Human qPCR Primer Pair (NM_005553)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3846
Forward Sequence:	TGTTCTGCTCTAGCTGTGGCA
Reverse Sequence:	AGGATGACCCACAGCCTGAAGA
ACCN:	<u>NM_005553</u> , <u>NM_005553.1</u> , <u>NM_005553.2</u> , <u>NM_005553.3</u> , <u>BC101744</u> , <u>BC101744.1</u> , <u>BC069531</u> , <u>NM_005553.4</u>
UniProt ID:	<u>P26371</u>
Synonyms:	KRN1; KRTAP5-1; KRTAP5.9
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)