

## Product datasheet for **HP208671**

### **KRT6A Human qPCR Primer Pair (NM\_005554)**

#### **Product data:**

<b>Product Type:</b>	qPCR Primer Pairs
<b>Gene ID:</b>	3853
<b>Forward Sequence:</b>	GAGGAGATTGCTCAGAGAAGCC
<b>Reverse Sequence:</b>	CAATCTCCTGCTTGGTGTTCG
<b>ACCN:</b>	<u><a href="#">NM_005554</a></u> , <u><a href="#">NM_005554.1</a></u> , <u><a href="#">NM_005554.2</a></u> , <u><a href="#">NM_005554.3</a></u> , <u><a href="#">BC014152</a></u> , <u><a href="#">BC014152.2</a></u> , <u><a href="#">BC008807</a></u> , <u><a href="#">BC069269</a></u> , <u><a href="#">BC125058</a></u> , <u><a href="#">BC139753</a></u> , <u><a href="#">BG059253</a></u> , <u><a href="#">BT006899</a></u> , <u><a href="#">BT007605</a></u> , <u><a href="#">NM_005554.4</a></u>
<b>UniProt ID:</b>	<u><a href="#">P02538</a></u>
<b>Synonyms:</b>	CK-6C; CK-6E; CK6A; CK6C; CK6D; K6A; K6C; K6D; KRT6C; KRT6D; PC3
<b>Components:</b>	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
<b>Quality Control:</b>	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
<b>Storage:</b>	Store at -20°C.
<b>Stability:</b>	The primer mix is stable for one year from date of shipping.



[View online »](#)