

Product datasheet for **HP208673**

Cytokeratin 7 (KRT7) Human qPCR Primer Pair (NM_005556)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3855
Forward Sequence:	TGTGGATGCTGCCTACATGAGC
Reverse Sequence:	AGCACCACAGATGTGTCGGAGA
ACCN:	NM_005556 , NM_005556.1 , NM_005556.2 , NM_005556.3 , BC002700 , BC002700.2 , BC009521 , BC042076 , BC107082 , BC107083 , BE740515 , BE744321 , BU527162
UniProt ID:	P08729
Synonyms:	CK7; K2C7; K7; SCL
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)