

Product datasheet for **HP208682**

SLP76 (LCP2) Human qPCR Primer Pair (NM_005565)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3937
Forward Sequence:	GGAAGAAGCCACCTGTGCCAAA
Reverse Sequence:	GTCATAGGAAGTAGTGCTGGC
ACCN:	NM_005565 , NM_005565.2 , NM_005565.3 , NM_005565.4 , BC016618 , BC016618.1 , BG911397 , BT007273 , BY995095 , NM_005565.5
UniProt ID:	Q13094
Synonyms:	IMD81; SLP-76; SLP76
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)