

Product datasheet for **HP208683**

LDHA Human qPCR Primer Pair (NM_005566)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3939
Forward Sequence:	GGATCTCCAACATGGCAGCCTT
Reverse Sequence:	AGACGGCTTTCTCCCTCTTGCT
ACCN:	NM_005566 , NM_005566.1 , NM_005566.2 , NM_005566.3 , BC067223 , BC067223.1 , BG424022 , BG760757 , BM995301 , BU507758 , NM_005566.4
UniProt ID:	P00338
Synonyms:	GSD11; HEL-S-133P; LDHM; PIG19
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)