

## Product datasheet for **HP208684**

### LGALS3BP Human qPCR Primer Pair (NM\_005567)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3959
Forward Sequence:	GGCAGCAATGTCACCATGAGTG
Reverse Sequence:	GTGGAAGCACTTGACTGACGAC
Accession No:	<u><a href="#">NM_005567</a></u> , <u><a href="#">BC015761</a></u> , <u><a href="#">BC002403</a></u> , <u><a href="#">BC002998</a></u> , <u><a href="#">BQ649131</a></u>
Synonyms:	90K; BTBD17B; CyCAP; gp90; M2BP; MAC-2-BP; TANGO10B
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix is stable for one year from date of shipping. Store at -20°C.



[View online »](#)