

Product datasheet for **HP208705**

Nuclear Factor 1 (NFIA) Human qPCR Primer Pair (NM_005595)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4774
Forward Sequence:	GTGGAGGATGAAATGGACAGTCC
Reverse Sequence:	CTGCTGAAACCAGACTTCTCCG
ACCN:	<u>NM_005595</u> , <u>NM_005595.1</u> , <u>NM_005595.2</u> , <u>NM_005595.3</u> , <u>NM_005595.4</u> , <u>BC022264</u> , <u>BX537823</u> , <u>NM_005595.5</u>
UniProt ID:	<u>Q12857</u>
Synonyms:	BRMUTD; CTF; NF-I/A; NF1-A; NFI-A; NFI-L
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)