

Product datasheet for **HP208713**

POU3F2 Human qPCR Primer Pair (NM_005604)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5454
Forward Sequence:	GTGTTCTCGCAGACCACCATCT
Reverse Sequence:	GCTGCGATCTTGCTATGCTCG
ACCN:	<u>NM_005604</u> , <u>NM_005604.1</u> , <u>NM_005604.3</u> , <u>BC051699</u> , <u>BC051699.1</u> , <u>BQ423615</u> , <u>BU183021</u> , <u>NM_005604.4</u>
UniProt ID:	<u>P20265</u>
Synonyms:	brn-2; BRN2; N-Oct3; oct-7; OCT7; OTF-7; OTF7; POUF3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)