

## Product datasheet for **HP208718**

### **RbAp48 (RBBP4) Human qPCR Primer Pair (NM\_005610)**

#### **Product data:**

|                          |  |
|--------------------------|--|
| <b>Product Type:</b>     | qPCR Primer Pairs  |
| <b>Gene ID:</b>          | 5928   |
| <b>Forward Sequence:</b> | AGACTTGCGTCTCCGTGGACAT   |
| <b>Reverse Sequence:</b> | CCTCCTTTGGAACGGCACTGAT   |
| <b>ACCN:</b>             | <u><a href="#">NM_005610</a></u> , <u><a href="#">NM_005610.1</a></u> , <u><a href="#">NM_005610.2</a></u> , <u><a href="#">BC053904</a></u> , <u><a href="#">BC053904.1</a></u> , <u><a href="#">BC075836</a></u> , <u><a href="#">BC003092</a></u> , <u><a href="#">BC015123</a></u> , <u><a href="#">BT007309</a></u>                     |
| <b>UniProt ID:</b>       | <u><a href="#">Q09028</a></u>  |
| <b>Synonyms:</b>         | lin-53; NURF55; RBAP48   |
| <b>Components:</b>       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.   |
| <b>Quality Control:</b>  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| <b>Storage:</b>          | Store at -20°C.  |
| <b>Stability:</b>        | The primer mix is stable for one year from date of shipping.   |



[View online »](#)