

Product datasheet for **HP208727**

ACSM3 Human qPCR Primer Pair (NM_005622)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6296
Forward Sequence:	CTGGGCAAAGTCTGCATGGAGT
Reverse Sequence:	AGTTGGTGCTGAACAGAAGACTG
ACCN:	<u>NM_005622</u> , <u>NM_005622.1</u> , <u>NM_005622.2</u> , <u>NM_005622.3</u> , <u>BC002790</u> , <u>BC002790.2</u> , <u>BI770558</u> , <u>BI826935</u> , <u>BM828310</u> , <u>BU170035</u> , <u>NM_005622.4</u>
UniProt ID:	<u>Q53FZ2</u>
Synonyms:	SA; SAH
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)