

Product datasheet for **HP208732**

SLC1A5 Human qPCR Primer Pair (NM_005628)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6510
Forward Sequence:	TCCTCTTCACCCGCAAAAACCC
Reverse Sequence:	CCACGCCATTATTCTCCTCCAC
ACCN:	<u>NM_005628</u> , <u>NM_005628.1</u> , <u>NM_005628.2</u> , <u>BC000062</u> , <u>BC000062.2</u> , <u>NM_005628.3</u>
UniProt ID:	<u>Q15758</u>
Synonyms:	AAAT; ASCT2; ATBO; M7V1; M7VS1; R16; RDRC
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)