

Product datasheet for **HP208744**

TAF6 Human qPCR Primer Pair (NM_005641)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6878
Forward Sequence:	CCTGATGCGTATGGTGAAAGCG
Reverse Sequence:	ACATCTGGTCGCAGGCACAACT
ACCN:	<u>NM_005641</u> , <u>NM_005641.1</u> , <u>NM_005641.2</u> , <u>NM_005641.3</u> , <u>BC018115</u> , <u>BC018115.1</u> , <u>BC127819</u> , <u>NM_005641.4</u>
UniProt ID:	<u>P49848</u>
Synonyms:	ALYUS; MGC:8964; TAF(II)70; TAF(II)80; TAF2E; TAFII-70; TAFII-80; TAFII70; TAFII80; TAFII85
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)