

## Product datasheet for **HP208746**

### TAF11 Human qPCR Primer Pair (NM\_005643)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6882
Forward Sequence:	TCCATCACTGGCACCTCTGTGT
Reverse Sequence:	GGTTGTAGTGGTGGCATTCTCC
ACCN:	<a href="#">NM_005643</a> , <a href="#">NM_005643.2</a> , <a href="#">BC021972</a> , <a href="#">BC021972.2</a> , <a href="#">BT019404</a> , <a href="#">BT019405</a> , <a href="#">BX538195</a> , <a href="#">BX647568</a> , <a href="#">NM_005643.4</a>
UniProt ID:	<a href="#">Q15544</a>
Synonyms:	MGC:15243; PRO2134; TAF2I; TAFII28
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)