

Product datasheet for **HP208748**

TAF13 Human qPCR Primer Pair (NM_005645)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6884
Forward Sequence:	GATGACCAGAATCCTTATACTGAG
Reverse Sequence:	ACTTGACTCGACCTTGTCTTCC
ACCN:	<u>BC017821</u> , <u>NM_005645</u> , <u>NM_005645.1</u> , <u>NM_005645.2</u> , <u>NM_005645.3</u> , <u>BC121180</u> , <u>BC121181</u> , <u>BG675152</u> , <u>NM_005645.4</u>
UniProt ID:	<u>Q15543</u>
Synonyms:	MRT60; TAF(II)18; TAF2K; TAFII-18; TAFII18
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)