

Product datasheet for **HP208764**

EVI5 Human qPCR Primer Pair (NM_005665)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7813
Forward Sequence:	AGGCTTCAGGAAGAACTCATTGC
Reverse Sequence:	CCAGTAGTACGAGCTAAGTGCC
ACCN:	<u>NM_005665</u> , <u>NM_005665.1</u> , <u>NM_005665.2</u> , <u>NM_005665.3</u> , <u>NM_005665.4</u> , <u>NM_005665.5</u> , <u>BC141544</u> , <u>BC156793</u> , <u>BP282704</u> , <u>NM_005665.6</u>
UniProt ID:	<u>O60447</u>
Synonyms:	EVI-5; NB4S
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)