

Product datasheet for **HP208783**

FARSLB (FARSB) Human qPCR Primer Pair (NM_005687)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10056
Forward Sequence:	GCTTACCGAAAGGAGATGGTGAG
Reverse Sequence:	CTGATTCCCATCACCTATGACTTC
ACCN:	<u>NM_005687</u> , <u>NM_005687.1</u> , <u>NM_005687.2</u> , <u>NM_005687.3</u> , <u>NM_005687.4</u> , <u>BC017783</u> , <u>BC006502</u> , <u>BC008434</u> , <u>BC011261</u> , <u>BC022273</u> , <u>BG721840</u> , <u>NM_005687.5</u>
UniProt ID:	<u>Q9NSD9</u>
Synonyms:	FARSLB; FRSB; HSPC173; NEDBLLA; PheHB; PheRS; RILDBC; RILDBC1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)