

Product datasheet for **HP208788**

LXR alpha (NR1H3) Human qPCR Primer Pair (NM_005693)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10062
Forward Sequence:	TGGACACCTACATGCGTCGCAA
Reverse Sequence:	CAAGGATGTGGCATGAGCCTGT
ACCN:	<u>NM_005693</u> , <u>NM_005693.1</u> , <u>NM_005693.3</u> , <u>BC008819</u> , <u>BC041172</u> , <u>BM835230</u> , <u>BP300071</u> , <u>BT019411</u> , <u>NM_005693.4</u>
UniProt ID:	<u>Q13133</u>
Synonyms:	LXR-a; LXRA; RLD-1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)