

Product datasheet for **HP208793**

GTP binding protein era homolog (ERAL1) Human qPCR Primer Pair (NM_005702)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	26284
Forward Sequence:	CTGGACCACTTCCTCGGATTCT
Reverse Sequence:	TCAGGGTGATGGACCAAGAGGA
ACCN:	NM_005702 , NM_005702.1 , NM_005702.2 , NM_005702.3 , BC019094 , BC019094.2 , BI758329 , BM683046 , BX384612 , BX537542 , NM_005702.4
UniProt ID:	O75616
Synonyms:	CEGA; ERA; ERA-W; ERAL1A; ERAL1B; H-ERA; HERA-A; HERA-B; PRLTS6
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)