

Product datasheet for **HP208808**

ARPC1B Human qPCR Primer Pair (NM_005720)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10095
Forward Sequence:	CTCCTGTGACTTCAAGTGTCGG
Reverse Sequence:	GGCTGAGAAACAGACGCCATGT
ACCN:	<u>NM_005720</u> , <u>NM_005720.1</u> , <u>NM_005720.2</u> , <u>NM_005720.3</u> , <u>BC002988</u> , <u>BC002988.1</u> , <u>BC002562</u> , <u>BC007555</u> , <u>NM_005720.4</u>
UniProt ID:	<u>O15143</u>
Synonyms:	ARC41; IMD71; p40-ARC; p41-ARC; PLTEID
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)