

Product datasheet for **HP208820**

ACTR1B Human qPCR Primer Pair (NM_005735)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10120
Forward Sequence:	GGTTGACTTCATACCTCGGCT
Reverse Sequence:	TCTGGCAACGTGTACTGCACCT
ACCN:	NM_005735 , NM_005735.1 , NM_005735.2 , NM_005735.3 , BC004374 , BC006372 , BC010090 , BC025274 , NM_005735.4
UniProt ID:	P42025
Synonyms:	ARP1B; CTRN2; PC3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)