

Product datasheet for **HP208833**

AKAP9 Human qPCR Primer Pair (NM_005751)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10142
Forward Sequence:	GGCGTCATTGATGGCTATGCAG
Reverse Sequence:	GCTGTTGCTCTGCCTCCAATTC
ACCN:	<u>NM_005751</u> , <u>NM_005751.1</u> , <u>NM_005751.2</u> , <u>NM_005751.3</u> , <u>NM_005751.4</u> , <u>BC015533</u> , <u>BC027455</u> , <u>BC040932</u> , <u>BC062448</u> , <u>NM_005751.5</u>
UniProt ID:	<u>Q99996</u>
Synonyms:	AKAP-9; AKAP350; AKAP450; CG-NAP; HYPERION; LQT11; MU-RMS-40.16A; PPP1R45; PRKA9; YOTIAO
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)