

Product datasheet for **HP208838**

ABI2 Human qPCR Primer Pair (NM_005759)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10152
Forward Sequence:	CAGCCAACCTTGAACGACCAGT
Reverse Sequence:	CTGAGTTGGAGGTGTTGTACGC
ACCN:	NM_005759 , NM_005759.1 , NM_005759.2 , NM_005759.3 , NM_005759.4 , NM_005759.5 , BC001439 , BC001439.2 , BC033944 , BT009920 , NM_005759.7
UniProt ID:	Q9NYB9
Synonyms:	ABI-2; ABI2B; AbIBP3; AIP-1; AIP1; argBP1; argBPIA; argBPIB; SSH3BP2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)