

Product datasheet for **HP208849**

RNA terminal phosphate cyclase like 1 (RCL1) Human qPCR Primer Pair (NM_005772)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10171
Forward Sequence:	CCTGTGAGGAAGGCTTGAAGC
Reverse Sequence:	GTTCGCCATCTGAGGTGACACA
ACCN:	NM_005772 , NM_005772.1 , NM_005772.2 , NM_005772.3 , NM_005772.4 , BC001025 , BC001025.2 , BI859171 , NM_005772.5
UniProt ID:	Q9Y2P8
Synonyms:	RNAC; RPCL1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)